

Food Rheology Applications

Fritz Soergel, Annika Hodapp, Philipp Beutler, Klaus Oldoerp, Fabian Meyer
Application Specialists, Material Characterization, Karlsruhe, Germany

17 May 2022

 The world leader in serving science



Contents

Rheology, Texture, Tribology ... (possible with one rheometer)

Rheology

- **Bread spreads, cream cheese, mayonnaise** (spreadability) Vane
- **Chocolate spreads, mayonnaise** (viscosity curve, stress plateau) Plate/cone
- **Chocolate melt** (ICA standard method 46 International Confectionery Association) Concentric cylinder
- Fresh fruit **juices**, tomato **sauce** with tomato & onion pieces Plate/plate (large gap)
- **Starch** (temperature induced transition crystallin/amorphous) Vane, RheoScope

Texture

- **Marshmallows** (chewing simulation) Plate
- Loaf **bread** slices (TPA and firmness) Plate
- **Bread roll**, apple, **margarine** (penetration, indenter test) Indenter
- **Chocolate, Cookies** (breaking test, crunchiness) Indenter + Lower

Tribology

- **Butter Milk Products** Ball-on-3-pins

Rheology & Texture

- **Cheese** (milk-based and analogue vegan) Serrated plate/p late, vane



Bread spread, cream cheese, mayonnaise (spreadability)



CR-Rotation vane rotor testing in original container holder

Chocholate spread	τ_0 in Pa	t in s
A	364	117
B	722	73

Tab. 2: Comparison of FL 22 vane rotor yield stress measurements with 0.05 rpm for chocolate hazelnut spread products




Fig. 5: FL 22 vane rotor yield stress measurements with 0.05 rpm for chocolate hazelnut spread products A and B

Raspberry jam	τ_0 in Pa	t in s
With seeds	590	100
Sieved	967	104
Sieved (measurement 2)	977	104

Tab. 1: Comparison of FL 22 vane rotor yield stress measurements with 0.05 rpm for a sieved raspberry jam and one with seeds




Fig. 4: Comparison of FL 22 vane rotor yield stress measurements with 0.05 rpm for a sieved raspberry jam (blue curve) and raspberry jam with seeds (orange curve, inserted image)

Fat content in %	τ_0 in Pa at 0.02 rpm
22.5	95
10.5	126
5.2	144

mayonnaises with three different fat contents




Fig. 6: Vane rotor yield stress measurements with 0.02 rpm for mayonnaises with three different fat contents

APPLICATION NOTE

No. V238

Spreadability of cream cheese - Influence of temperature and fat content

APPLICATION NOTE

No. V271

Investigation the effect of fat content on the yield stress of mayonnaise measured with HAAKE Viscotester iQ Rheometer and vane rotor

APPLICATION NOTE

No. V272

Yield stress of jam, chocolate spread and peanut butter measured with HAAKE Viscotester iQ Rheometer and vane rotors




- Thermo Scientific Rheology Application Note V238 Applied Food Rheology Using Fast Speed Control and Axial Measurements
- Thermo Scientific Rheology Application Note V268 Spreadability of Cream Cheese - Influence of Temperature and Fat Content
- Thermo Scientific Rheology Application Note V271 Investigation of the effect of fat content on the yield stress of mayonnaise
- Thermo Scientific Rheology Application Note V272 Yield stress of jam, chocolate spread and peanut butter

Chocolate spread, mayonnaise (viscosity curve, stress plateau)


CR-Rotation equilibrium steps  (are more exact than ramps ), air bearing rheometers

- **Plate/cone** measuring geometry **temperature-controlled**
- Viscosity decreases by 7 or 3 orders of magnitude, resp.

Viscosity curve of a chocolate spread.
MARS 60 with **plate/cone** C35/1° at 35 °C



Mayonnaise type	Steady-state step-test stress plateau
Regular	116 Pa
Low fat	89 Pa
40% yogurt	119 Pa



thermoscientific

APPLICATION NOTE

Rheological and textural properties of various food formulations analyzed with a modular rheometer setup

ANIS163

Authors
Fabian Meyer and Klaus Oldörp
Thermo Fisher Scientific, Karlsruhe, Germany

Key words
Rheology, food, viscosity, yield stress, texture analysis, tribology

Introduction
Since foodstuff comes in such a vast variety of structures and textures, the range of rheological methods used to characterize its mechanical properties is even wider.

Rheological and texture properties play an important role during the entire life cycle of liquid or solid food formulations. Starting with simple single point viscosity measurements in the original containers for batch release in production, over the determination of classical rheological parameters like shear viscosity or yield point for quality control purposes, the mechanical testing of foods reaches a certain level of complexity with comprehensive rheological investigations for the development of new formulations in the research and development department.

While some methods rely on classic rheometer geometries like parallel plates, cone & plate or coaxial cylinders, some other methods try to emulate a certain application by utilizing special rotors and/or measuring fixtures. One such application is studying the texture of food products, which has to match the consumer's expectations. With specially designed probes a rotational rheometer can test various important food properties such as softness, stickiness or spreadability, and can even be utilized for performing tribological tests.

In this application selected rheometer accessories and physical measurements of various food products will be reviewed using a modular rheometer designed for product development and quality control purposes. This includes measuring "classic" rheological properties such as viscosity and yield stress, as well as customized setups for a comprehensive mechanical investigation of food formulations.

Materials and methods
Different commercially available food products were examined using the Thermo Scientific™ HAAKE™ MARS™ IQ Rheometer with a mechanical bearing (Figure 1). The tests to be performed included viscosity and yield stress measurements, axial bending, breaking and squeezing tests as




Figure 1: HAAKE MARS IQ Rheometer with Peltier plate temperature control for use with parallel-plate and cone-and-plate geometries.

ThermoFisher
SCIENTIFIC



Chocolate (ICA standard method 46 International Confectionery Association)



CR-Rotation ramp, thixotropy loop area, curve fits and interpolation for yield stress determination

- **Concentric Cylinder** measuring geometry @ 40°C
- **Push Button** for sample conditioning, measurement, evaluation & documentation

	Milk Chocolate	Dark Chocolate
τ_0 Casson / Pa	8.9	2.1
τ_0 Windhab / Pa	14.7	4.0
τ_0 Servais et al. / Pa	30.0	10.4



thermo scientific

APPLICATION NOTE

Flow behaviour of chocolate melts-working according to ICA standards

Author: Klaus Oldörp

Keywords
Chocolate, QC, ICA standards, viscosity, yield stress

Introduction
The flow behaviour of molten chocolate is a crucial parameter for many reasons. During production the transport, filling, dipping, coating or dosing steps depend on a well defined viscosity and yield stress. Likewise, the properties of the final chocolate like the look of its surface or its mouth feeling are directly related to the chocolate's viscous behaviour.

Testing the viscosity is therefore one of the standard quality control (QC) test methods for any company producing chocolate or using chocolate for their own production of e.g., chocolate-coated cookies.

To make viscosity testing in QC easier and more reliable, the Thermo Scientific™ HAAKE™ Viscotester™ IQ Rheometer (Figure 1) has been developed. This viscometer includes features especially designed for QC applications. For example, due to its improved sensitivity it is possible to use smaller measuring geometries, which reduces sample volume, time for temperature equilibration and cleaning effort. Also, with improved sensitivity smaller shear rates are accessible, which improves the reliability of yield stress calculations [1] with extrapolation methods like the Casson model.

Preparations
Two chocolate samples, a milk chocolate and a dark chocolate, have been prepared according to ICA method 46 [2] by putting chocolate pieces into glass containers, sealing the containers and leaving them in an oven at 52 °C for between 45 and 60 minutes. Meanwhile the cup and bob of the measuring geometry are preheated to 40 °C the HAAKE Viscotester IQ Rheometer in the Pelletier temperature control unit.





ThermoFisher SCIENTIFIC

Fresh fruit juices, tomato sauce with tomato & onion pieces

CR-Rotation viscosity curve: large gap **plate/plate** meas. geom. TMP60 DC (10 mm rim) + P60 rotor

CS-Rotation ramp for yield stress measurement and stability prediction

- Due to large particles (several mm) absolute measurements are not possible
- Fresh pineapple juice with fruit particles (4 mm gap)
- Tomato sauce with larger tomato & onion pieces (10 mm gap)



Starch (temperature sweeps)

CR-Rotation with vane in temperature-controlled cup

CD-Oscillation plate/plate temperature sweeps with polarization light microscope

- **Temperature** induced transition crystallin/amorphous



Fig. 2: Propeller shaped 2-blades vane rotor FL26-2b with fitted shaft adapter (left) and 4-blades vane rotors FL22 and FL16 with shaft adapter (right)




Fig. 13: Wheat flour (type 405), two grinding grades, vane rotor FL26-2b; viscosity in arbitrary units




Fig. 15: Crushed rice and maize, same temperature profile as with wheat flour followed by a constant temperature of 65 °C, FL22 vane rotor; viscosity in arbitrary units



Fig. 14: Maize grain (front)






Figure 5: Native potato starch (5 % in water): Photos show the starch crystals in the beginning, the swollen crystals at peak viscosity and the inhomogeneous solution after cooling down.




Thermo Scientific Rheology Application Note V228 RheoOptical Tests during the cooking of starch and the crystallization of fats

Thermo Scientific Rheology Application Note V255 Rheological investigation of starch swelling, gelatinization and retrogradation

Milk, yoghurt, cheese, fat, dough (gelation, baking, crystallization)

CD-Oscillation time sweeps or temperature sweeps, **plate/plate** measuring geometry

- Oscillation
 - does not affect the development of the structure
 - Yields valuable information about the structure (e.g. strength)
- Gelling of protein (see graph)
- Baking of dough (e.g. flour batter)
- Crystallization of fat (cooling run)
- Solid elastic samples (e.g. cheese)



Marshmallows (chewing simulation)

Texture analysis with a plate

- **AXIAL** h-controlled lift movement




Fig. 1: Samples M (left) and S (right) before measurement.





Fig. 2: Samples M (left) and S (right) during measurement.

Repeated Compression of a Marshmallow



A Rheometer with Bite – Marshmallows in the Rheometer

Klaus Oldörp, Thermo Fisher Scientific, Material Characterization, Karlsruhe, Germany

Application Notes V243

Introduction

We use many of our senses to experience our food. For a food product to be popular and therefore successful, it must meet many requirements: it has to look appealing, smell good, have a pleasant texture in the mouth - or "mouth feel" - and of course taste good. As soon as a food product is perceived as being "unpleasant" during even one of these subjective "tests," it can quickly lead to consumer rejection.

When regarding natural products like bread, meat or cheese, certain fluctuating or individual "bad" qualities are sometimes accepted, such as the smell of certain cheeses. However for industrially-produced foodstuffs, and particularly when dealing with treats such as sweets, all of the above criteria have to be met to gain consumer acceptance. No matter how good something looks or tastes, the "mouth feeling," i.e. the texture in the mouth, has to be right, otherwise the pudding will remain on the shelf in the refrigerated section and the cookies will never be purchased.

Marshmallows are another example of food where not only the flavour and sweetness, but also the "bite", is critical for enjoyment. Using a rheometer or a testing machine such as a texture analyzer, this consistency or texture can be described using objective parameters and thus making quality control or targeted improvement possible. A modern rheometer allows various options for characterising marshmallows and similar products. It also allows the viscosities of the starting materials to be measured, and the finished product can be characterized in an oscillatory test.

By combining a sensitive normal force sensor with a high-precision lift drive, the Thermo Scientific™ HAAKE™ MARS™ provides the additional option of stressing samples axially, i.e. by exerting a vertical force on them from above at a maximum of 50 N (which corresponds to a weight of 5 kg) by either pressing on or pulling the marshmallow sample. In this type of situation, the Thermo Scientific HAAKE MARS measures the axial force and position of the measurement geometry precisely while the sample is being squeezed, or the penetration of a probe into the sample is tracked.




Fig. 1: Samples M (left) and S (right) before measurement.




Fig. 2: Samples M (left) and S (right) during measurement.

Measurements and Results

Two products were compared during the experiment. Marshmallows (M) of about 30 mm in height and a sugar foam product (S) of about 10 mm in height. In order to compare the absolute values discs were cut from product S matching the diameter of the marshmallows and were stacked three high (Fig. 1). A 35 mm-diameter plate as used for rotational and oscillation measurements was run down to a height of 7 mm (Fig. 2) at a speed of 1 mm/s, and then run back up to 30 mm at 1 mm/s. This process was repeated on the samples five times in a row to simulate chewing in the mouth.

When compressing first sample M, the force increased to about 45 N before dropping again. During the following cycles the force only reached about 40 N but remained nearly constant throughout each cycle. A force of 50 N is required to compress sample S; the sample is thus felt to be somewhat harder. The maximum applied force decreases throughout subsequent cycles, dropping from 44 N during the second cycle to 40 N during the fifth cycle.


Loaf bread slices (TPA and firmness)

Texture analysis with a plate

- **AXIAL** h-controlled lift movement


Texture profile analysis (TPA)

Two-bite test (2 compression cycles of a cylindrical sample by 4 mm (30 %) with 0.4 mm/s)



Firmness test based on AACC 74-09

Compression ε of two slices of loaf bread with 100 mm/s up to 40 % of the initial gap height and evaluation of compression force value (CFV) at $\varepsilon = 25$ %



Bread roll, apple, margarine (penetration, indenter test)

Texture analysis with an indenter

- **AXIAL** h-controlled lift movement

Maximum force	t in s	F in N
Upper half	8.2	6.0
Lower half	5.8	4.6



Crust : $F_{\text{Upper half}} > F_{\text{Lower half}}$

Crumb : Compression of the inner part


Indenter diameter 6 mm

Lift speed 1.2 mm/s

Penetration test on margarine repeated 8 times.
First the probe moves into the sample,
Then the relaxation of the normal force is measured.
Indenter diameter 5 mm



Upper half




Lower half

Chocolate (breaking test)

Texture Analysis – Chocolate experts can judge the quality of a chocolate bar by breaking off a piece

- **AXIAL** h-controlled lift movement with 1.3 mm/min, $T = 23^{\circ}\text{C}$
- **Indenter** diameter 5 mm
- Force for breaking correlates with sound of empirical breaking test



thermo scientific

APPLICATION NOTE

Breaking strength of chocolate

Authors: Sabine Dohmen, Klaus Oldörp, Hans-Michael Petri

Chocolate appeals to all of our senses. We see the colour and the glossy surface of the bar, we smell the elusive flavours, we weigh it in our hands and feel how it melts, we hear it break and sense its resistance, and savour the taste. With so much "sensuality", it is easy to forget how quickly the pleasure can be destroyed, if even just one of these properties is not as we expect it or are accustomed to.

Even unconsciously, the breaking behaviour of chocolate plays an important role in influencing the consumer's impression. Chocolate experts can even evaluate the quality of a sample by breaking off a piece.

Rheology plays an important role in several steps of chocolate production. The liquid chocolate formulation and the fats used can be characterized by their viscosities, yield stresses and solidification behaviours. These parameters are important for quality control and processing, and can be determined using rotational or oscillation measurements.



Sensory properties as experienced during the melting process, or "mouth feeling", can be described with the viscosity curves and the yield stress. However, these rheological parameters contain no information about the breaking strength of the final chocolate bar – and therefore a new method and new equipment are necessary to assess the resistance to break.

The Thermo Scientific™ HAAKE™ MARS™ Rheometer, manufactured by Thermo Fisher, features a highly sensitive normal force sensor and a very precise lift motor which allows the customer to apply controlled axial forces to the sample, pushing or pulling it, and to analyze its axial deformation.

For example, with a new measuring geometry (Figure 3), chocolate bars can be positioned on the rheometer and submitted to an increasing axial force until they break. This measuring geometry [1] consists of two parallel support bars which can be mounted onto a base plate in a variable



Figure 1: The HAAKE MARS Rheometer.

Cookies (breaking test, crunchiness)

Texture Analysis

- **AXIAL** h-controlled lift movement, **indenter**
- Fresh cookies (reference)
- Aged cookies (exposed to air/humidity) more elastic
 - Good Reproducibility
 - 22 % higher force F_{nmax}
 - 8 % more bending Δh_{max}





Table 1: Amount of bending and force at curve maximum of fresh and aged biscuits.

No.	Fresh biscuits		Biscuits aged for 14 day	
	Δh / mm	F_n / N	Δh / mm	F_n / N
1	0.532	5.384	0.633	6.395
2	0.563	6.395	0.596	7.927
3	0.596	5.943	0.584	7.545
4	0.501	5.332	0.530	7.289
5	0.490	6.684	0.562	7.065
Mean value	0.537	5.948	0.581	7.244
StdDev	0.039	0.536	0.034	0.512
rel. StdDev	7.31 %	9.02 %	5.90 %	7.07 %



thermoscientific

APPLICATION NOTE

Crunchiness of biscuits – Texture analysis on a rheometer

No. V293

Author
Klaus Oldörp,
Thermo Fisher Scientific, Karlsruhe, Germany

Key words
Cookies, biscuits, crunchiness, texture analysis, bending and breaking, rheometer

Introduction
Biscuits or cookies, no matter what they are called, they need to have the right taste and the right bite. While the first is a rather subjective matter and can lead to lengthy discussions, the latter can be tested in a very objective way with a rheometer.

Material and methods
For testing biscuits, a Thermo Scientific™ HAAKE™ MARS™ iQ Rheometer (Figure 1) was equipped with the 3-point bending-accessory [1], using the 8 mm plate as the probe (Figure 2). Other probes can easily be adapted to the measuring head using the universal adapters U1 with a 6 mm bore or U2 with a 4 mm bore.

Apart from the classic measuring parameters shear rate, deformation and stress, a modern rheometer comes equipped with an additional sensor to detect normal force. Thus, combined with its precise lift a rheometer can do texture analysis as well.

The measuring protocol starts with zeroing the normal force and lowering the probe to a position still high enough to allow a quick placement of a fresh biscuit onto the lower supports of the 3-point-bending accessory. For the tests discussed here, a position about 2 mm above the biscuits upper surface has been used.


The lower supports of the 3-point-bending-tool were positioned 5 cm apart from each other to support the biscuit approximately 0.8 cm from its edges. The biscuit was placed onto the lower supports with its centre under the probe (Figure 2).

Figure 1: The Thermo Scientific HAAKE MARS iQ Rheometer.

ThermoFisher
SCIENTIFIC

Different Tribology measuring geometries available – here: **Ball-on-3 Pins** geometry

- 3D printed (rapid prototyping), machined version available
- **Rotor** with ceramic shaft, metal spring and holder for a ½" **ball** (here: glass)
- **Lower holder** for **3 Pins** (here: PDMS, made in a 3D printed mould)



Thermo Scientific Rheology Application Note (in prep.) Influence of fat content on the tribological properties of butter milk products

Shewan H.M., Pradal C, Stokes J.R.: Tribology and its growing use toward the study of food oral processing and sensory perception. J. Texture Stud. 2020-51 7–22

Cheese – one rheometer for rheology and texture testing

Comparison milk-based cheese vs. vegan cheese analogue

Serrated plates temperature controlled

- Normal force relaxation (3 N initially)
- CD-oscillation amplitude sweep
- CD-oscillation temperature sweep
- Texture serrated plate

Vane rotor

- Cutting/biting test with vane rotor blades



Figure 2: Upper and lower serrated plate (left) and measuring geometry with cheese sample (right).



Figure 3: 4-blade vane rotor used to simulate cutting or biting the product.




Figure 7: Results from the amplitude sweeps at 37 °C plotted over the deformation γ in %.




Figure 8: Development of viscosity and phase angle during heating up to 90 °C.




Figure 9: Results of penetration test of both products.



Rheology and texture of cheese and a vegan cheese analogue

Authors

Annika Hodapp and Klaus Oldorp, Thermo Fisher Scientific, Karlsruhe, Germany

Keywords

Cheese, vegan cheese analogue, food analogue, texture analysis, rheology, rheometer

Introduction

What is cheese? Cheese is fermented milk with the right taste and the right texture. While the taste is far beyond its capabilities, a rheometer can definitely help to analyse the texture. This becomes especially helpful when trying to create a vegan substitute for cheese, since the consumer expects a similar look and feel compared to a milk-based cheese.

The texture is an important part of a cheese's properties. Depending on the kind of cheese and its age this can range from creamy to solid. Rheological measurements can be used to quantify the viscoelastic properties of cheese and cheese formulations. On top of that, a modern rheometer is equipped with a normal force sensor and therefore offers capabilities beyond the classic rheological measurements. In combination with its automatic lift, a rheometer can also be used to run texture analysis tests. In the case of a cheese this could be used to quantify e.g. its bite or its cutting properties.

This report summarises the results of various rheological measurements and a penetration test on a milk-based cheese and a vegan cheese substitute.

Materials and methods

For the measurements described here, a sliced butter cheese and a sliced vegan alternative were used. Both products are commercially available and were purchased at a local grocery shop. The milk-based cheese consisted mainly of protein and fat (25 % and 24 % respectively). The main components of the vegan alternative were fat and starch (20 % and 17 % respectively) while its protein content was below 1 %.




Figure 1: The Thermo Scientific HAAS MARS iQ Rheometer.

The rheological tests were performed using a Thermo Scientific™ HAAS™ MARS™ iQ Air Rheometer (Figure 1) equipped with a Peltier temperature module for cone-plate- or parallel-plate geometries. To be able to adapt to the thickness of the cheese slices, a parallel-plate-geometry was used. To prevent slippage, serrated plates were used on both sides (Figure 2).

To get a general idea about the properties of both products, amplitude sweeps at 20 °C and 37 °C were performed over a strain range from 0.01 % and 100 % at a frequency of 1 Hz. The measurements at 20 °C were intended to show the sample properties at room temperature, where cheese is cut or a piece is bitten off, for example. At 37 °C the behaviour during chewing or swallowing was investigated.

Conclusion

Rheology, Texture, Tribology ... (possible with one rheometer)

Rheology

- | | |
|--|-------------------------|
| • Bread spreads, cream cheese, mayonnaise (spreadability) | Vane |
| • Chocolate spreads, mayonnaise (viscosity curve, stress plateau) | Plate/cone |
| • Chocolate melt (ICA standard method 46 International Confectionery Association) | Concentric cylinder |
| • Fresh fruit juices , tomato sauce with tomato & onion pieces | Plate/plate (large gap) |
| • Starch (temperature induced transition crystallin/amorphous) | Vane, RheoScope |

Texture

- | | |
|--|------------------|
| • Marshmallows (chewing simulation) | Plate |
| • Loaf bread slices (TPA and firmness) | Plate |
| • Bread roll , apple, margarine (penetration, indenter test) | Indenter |
| • Chocolate, Cookies (breaking test, crunchiness) | Indenter + Lower |

Tribology

- | | |
|-------------------------------|----------------|
| • Butter Milk Products | Ball-on-3-pins |
|-------------------------------|----------------|

Rheology & Texture

- | | |
|---|-----------------------------|
| • Cheese (milk-based and analogue vegan) | Serrated plate/p late, vane |
|---|-----------------------------|



Thermo Scientific Food Rheology Application Notes, White Paper and Compendium

V136 Measurement and evaluation on liquid **chocolate** according to the IOCCC Standard

V187 Rheological behavior of hydro**gels** dependent on concentration of the polymer and type of neutralizing agent

V218 Testing **instrument** accuracy with Newtonian standard fluids

V222 Breaking strength of **chocolate**

V223 HAAKE **RheoWin** software features for quality control and routine measurements

V228 Rheo-optical Tests during the **cooking** of **starch** and the **crystallization** of **fats**

V238 Applied food rheology using fast speed control and **axial** measurements

V243 A rheometer with bite – **marshmallows** in the rheometer

V248 Well prepared – **good results**

V255 Rheological investigation of **starch** swelling, gelatinization and retrogradation

V260 Product **stability** or shelf life – what rheology has to do with it

V265 Image acquisition with the **RheoScope** module at very high shear rates (stroboscope, contrast enhancing illumination)

V268 **Spreadability** of cream cheese – influence of temperature and fat content

V269 Flow behaviour of **chocolate** melts – working according to ICA Standards

V271 Investigation of the effect of fat content on the yield stress of **mayonnaise**

V272 Yield stress of **jam**, **chocolate spread** and **peanut butter**

V275 How much squeezing power is required to get a **paste** out of a tube?

V292 Rheological and **textural** properties of various food formulations analyzed with a modular rheometer setup

V293 Crunchiness of biscuits – **texture** analysis on a rheometer

V297 Investigating **cocoa butter** crystallization using simultaneous rheology and **Raman** spectroscopy (RheoRaman)

V298 Investigating heat-induced **gelation of whey protein** using simultaneous rheology and FTIR spectroscopy

V300 Rheology and **texture** of **cheese** and a vegan cheese analogue

White Paper WP04 An evaluation of **meat analog product** characteristics by the combination of microscopy and rheological tools.

Compendium C001.2 Guide to skilled **food rheology & extrusion**

Thank you

