The new Büchi pressure filters – filtration on a new level

- easy operation
- removable filter basket
- quick opening / closing system
- perfect to Büchi pressure reactors
- for every-day lab use



Pilot Plant and Reactor Systems

The new pressure filter by Büchi – simple and superior

Accessories

- Various filter sizes and filter types available
- With or without heating jacket, adjusted to your needs



Small version up to a volume of 5 litres.



Large version up to a volume of 20 litres.



Clamp screws guarantee easy opening.







Simply lift the filter basket for removal of the filter cake and cleaning. A time-saving and cost efficient solution.

Visit www.buchiglas.com for more information about:

- pressure reactors
- glass plants
- automation