



## Kjelflex K-360

### Broad application flexibility

Perform Kjeldahl nitrogen or protein determination and additional steam distillation in the most flexible way. The modular setup allows to adapt the Kjelflex K-360 at any time including a wide selection of titrators to guarantee smooth and easy operation.

#### Convenient

Automate your analysis in the most supportive way



#### Data security

Improve documentation and data security with multi-level authorization



#### Flexible

Cope with the broadest range of applications for steam distillation

## K-360: Your most important benefits



### Convenient

- Customizable analysis methods for multi-purpose use
- Modular concept to adapt to changing needs
- Modern instrument software enabling convenient routine use
- Connection to various external titrators for semi-automated processing



### Data security

- Storage of results on the instrument (titrant consumption)
- Data export on USB printer for documentation purposes
- No unauthorized use due to password protection



### Flexible

- Acid resistant pump for additional steam distillation applications
- Wide selection of titrators for colorimetric or potentiometric procedures
- Reduction of steam power for gentle distillation
- Choice of glass or plastic splash protector

## Complete your portfolio



SpeedDigester  
K-439 / K-425 / K-436  
IR digestion



Extraction Unit  
E-816 SOX & HE  
Dedicated fat extraction



NIRMaster™  
IP54 / Pro IP65  
FT-NIR spectroscopy



Recirculating Chiller  
F-108  
The efficient way of  
cooling

