















Overview: Supported inkjet printers (1)

Supported printers: <ul style="list-style-type: none"> MI 9018, MI 9020, MI 9028, MI 9029, MI 9030 MI 9232, MI 94xx 	MI 9018	MI 9029	MI 9410	MI 9450
	 www.markem-ijaje.com/	 www.markem-ijaje.com/	 www.markem-ijaje.com/	 www.markem-ijaje.com/
Performance	<ul style="list-style-type: none"> 4 Lines Alphanumeric wide choice of characters Graphics/Logo printing Side wall printing 	<ul style="list-style-type: none"> 4 Lines Alphanumeric wide choice of characters Graphics/Logo printing Side wall printing 	<ul style="list-style-type: none"> 5 Lines Alphanumeric, bar codes, 2D codes Graphics/Logo printing Side wall printing 	<ul style="list-style-type: none"> 5 Lines Alphanumeric, bar codes, 2D codes Graphics/Logo printing Side wall printing
Other features	<ul style="list-style-type: none"> Bench top Key pad USB port RS-232 	<ul style="list-style-type: none"> Bench top Key pad USB port RS-232 	<ul style="list-style-type: none"> Bench top Touch screen USB and SD port Ethernet 	<ul style="list-style-type: none"> Bench top Touch screen USB and SD port Ethernet
Dimensions (HxWxD)	595 x 307 x ??? mm	595 x 307 x ??? mm	555 x 360 x 357 mm	555 x 360 x 357 mm
Weight (kg)	18.5	18.5	25	25
List Price (EUR)	<5'000	<10'000	>10'000	>10'000
Maintenance	<ul style="list-style-type: none"> ink exchange in < 6 min 2 year non-stop 	<ul style="list-style-type: none"> ink exchange in < 6 min ink runs > 8000 hours automatic cleaning system 	<ul style="list-style-type: none"> Maintenance: every 18 month Head cleaning recommended every second week Head can be left attached on MEDIAJET 	<ul style="list-style-type: none"> Maintenance: every 18 month Head cleaning recommended every second week Head can be left attached on MEDIAJET
INTEGRA comments	Industrial interface must be ordered.	Industrial interface must be ordered.	Recommended MI printer , lower consumable and maintenance costs. Industrial interface must be ordered.	Industrial interface must be ordered. Higher line speed and storage capacity.

Overview: Supported inkjet printers (2)

Supported printers:	Linx 89xx (89xx: -00, -10, -20, -40)	Domino Ax150i	Domino Ax350i
<ul style="list-style-type: none"> • Linx 4900, Linx 5900, Linx 89xx (printhead MK7) • Linx CJ400 (printhead MK11) • Domino A320i, Domino A420i, Domino Ax150i, Domino Ax350i 	 www.linxglobal.com	 www.domino-printing.com	 www.domino-printing.com
Performance	<ul style="list-style-type: none"> ▪ 6 Lines ▪ Alphanumeric, barcodes, 2D codes ▪ Graphics/Logo 25 px ▪ Side wall / bottom printing 	<ul style="list-style-type: none"> ▪ 5 Lines ▪ Alphanumeric, bar codes, 2D codes ▪ Graphics/Logo printing ▪ Side wall printing 	<ul style="list-style-type: none"> ▪ 5 Lines ▪ Alphanumeric, bar codes, 2D codes ▪ Graphics/Logo printing ▪ Side wall printing
Other features	<ul style="list-style-type: none"> ▪ Benchtop ▪ Touch screen ▪ USB port ▪ Ethernet 	<ul style="list-style-type: none"> ▪ Bench top ▪ Touch screen ▪ USB port ▪ Ethernet, RS232, GPIO 	<ul style="list-style-type: none"> ▪ Bench top ▪ Touch screen ▪ USB port ▪ Ethernet, RS232, GPIO
Dimensions (HxWxD)	503 x 440 x 355 mm	373 x 400 x 393 mm	411 x 430 x 331 mm
Weight (kg)	24	16	22
List Price (EUR)	>10'000	8'300	>10'000
Maintenance	<ul style="list-style-type: none"> ▪ Maintenance: every 24 months ▪ Head cleaning every 3 months 	<ul style="list-style-type: none"> ▪ Low-maintenance ▪ Head can be left attached on MEDIAJET 	<ul style="list-style-type: none"> ▪ Low-maintenance ▪ Head can be left attached on MEDIAJET
INTEGRA comments	This reliable Linx printer replaces all previous models.	Recommended Domino printer Recommended, if bar codes or 2D codes needed.	Recommended, if bar codes or 2D codes needed.

Ordering information MEDIAJET inkjet printer connections kits

Supported printers: <ul style="list-style-type: none"> MI 9018, MI 9020, MI 9028, MI 9029, MI 9030 MI 9232, MI 9410, MI 9450 Linx 4900, Linx 5900, Linx 89xx (print head MK7) Linx CJ400 (print head MK11) Domino A320i, Domino A420i, Domino Ax150i, Domino Ax350i 	Markem-Imaje (MI) 9018-9030	Markem-Imaje (MI) 9232/94xx	Linx side wall printing (with print head MK7 and MK11)	Linx bottom printing (with print head MK7 and MK11)	Domino
<ul style="list-style-type: none"> Linx 4900, Linx 5900, Linx 89xx (print head MK7) Linx CJ400 (print head MK11) Domino A320i, Domino A420i, Domino Ax150i, Domino Ax350i 					
Part number					
Printer connection kit: with printer head fixing device and interface cable 	103080	103091	113845	113840	113060
Operating instructions	Operating instructions for connection kits are included and also available in the download center on www.integra-biosciences.com				
	103983	103980	103956	103981	103982
Cart (W x D x H (cm): 70 x 50 x 65-90) 	-	-	-	113841	-
INTEGRA comments	Compatible with: Ø 90 mm, Ø 60 mm or Ø 35 mm Petri dishes.	Compatible with: Ø 90 mm or Ø 60 mm Petri dishes.	Compatible with: Ø 90 mm or Ø 60 mm Petri dishes.	Compatible with: Ø 90 mm, Ø 60 mm or Ø 35 mm Petri dishes.	Compatible with: Ø 90 mm or Ø 60 mm Petri dishes.
Compatibility with MEDIAJET (manufactured)	After Dec 2005 / Feb 2007; (SN above: MJ : 10305000 / MJ vario : 11307000)	After Feb 2013; (SN above: MJ : 10313568 / MJ vario : 11313247)	After Feb 2013; (SN above: MJ : 10313568 / MJ vario : 11313247)	After Mai 2010; (SN above: MJ : 10310371 / MJ vario : 11310128)	After Feb 2013; (SN above: MJ : 10313568 / MJ vario : 11313247)